

Step One


MiniOne Gel Kit Packing Guide Crate #1

Video Tutorial



bit.ly/HowToPackMini

Place **repair tags** on broken micropipettes.



Micropipettes and waterbath placed at the **bottom** with MiniOne supply box on **top** of micropipettes.

Step Two

MiniOne Gel Kit Packing Guide Crate #1

Video Tutorial



bit.ly/HowToPackMini

Stack tip
boxes **closed**
and upright.



Leave enough
space **above** tip
boxes so that
lid flaps **close**
completely with
no obstruction.



Science Education Partnership
fredhutch.org/SEP (206) 667-4487 sep@fredhutch.org